

The safety of your patients and their families is job #1.



Introducing the MinneMold™ 30-Cavity Troche Mold for professionally compounded troches

The MinneMold is an innovative troche mold design that enables pharmacists to dispense troches in a child-deterrent, tamper evident, easy to use system.



Easy to Use...

- Squeeze side tabs simultaneously while sliding tray out to dispense
- Slide tray in fully to lock

Easy to Compound...

- Fill tray cavities and allow troche base to cool. Insert tray into sleeve until it locks
- Apply optional tamper evident label
- Patent pending child deterrent design
- Tamper evident feature (labels included)
- Suitable for light-sensitive preparations (meets USP light transmission requirements)
- Made from an FDA approved resin
- BPA-Free polymer
- 30 count quarter-scored cavities
- 1.05g of Polyethylene Glycol 1450 per cavity



Only available from Letco Medical

MinneMold™
Clearly The Best.



facebook.com/letcomedical



linkedin.com/company/letco-medical



LETCO
MEDICAL

QUALITY COMPOUNDING STARTS HERE

800.239.5288 • www.letcomedical.com